


snappax

PACKAGING COLLECTION



Fold-to-open card

Folding the semi-rigid base opens the line at its centre. Pressure on the pack directs the product through this opening.



Printed flexible surface

Product reservoir and brand surface.

Semi-rigid reverse

Carrier folded along its centre.

Opening line

Fold the card along this line.



Physical sample: reverse and front

Fold. Squeeze. Apply.

Folding opens the intended line in the semi-rigid base. Pressure on the flexible film then dispenses the product through this opening.



01

Fold

Fold the card along its central opening line.



02

Squeeze

Press the flexible surface to dispense the contents.



03

Apply

Direct the product onto the area where it will be used.

Your needs define the shape.

We produce custom moulds to create single-dose packaging in your desired shape. Dimensions, cavities, filling and print are developed around your project.

Start with an example.

The products shown are examples, not a fixed selection. With custom mould production, we design the shape, dimensions and cavity layout around your requirements.



Fold-to-open card packs

Single-dose card packs that dispense the product through a central opening line.



Thermoformed cups

Custom-moulded single or multi-compartment cups sealed with lidding film.



Break-tip ampoule packs

Film-formed packs for liquid products, with a narrow neck and an integral break-off tip.



Break-tip ampoule

Weekly skin and hair care

Two seven-day collections for skin serum and hair serum. Each unit is marked in order from Monday to Sunday.



Multi-compartment pack

Six-step skincare set

Cleanse, tone, serum, eye care, moisturize and prime. Each of the two print alternatives has two columns and six rows.



Break-tip ampoule

Fragrance and care samples

Three fragrance alternatives alongside face serum, hair serum and care oil in single-use sample packs.



Fold-to-open card

Honey, face cream and tomato sauce

Fold the card along its centre line, then squeeze the chamber. Honey, face cream or tomato sauce dispenses through the central opening.



Break-tip ampoule

Olive oil and balsamic vinegar

Olive oil and balsamic vinegar in transparent reservoirs that reveal the contents. Each three-unit strip contains individual portions with narrow necks and integrated tips.



Thermoformed cup

Masks, cream and cleansing gel

Face mask, hair mask, face cream and cleansing gel share the same cavity structure. Opaque printed surfaces are combined with clear bodies that reveal the fill.



Thermoformed cup

Mango and avocado purée

Mango and avocado purée presented in single-portion packs, with the fill contained in a moulded cavity.



Thermoformed cup

Eye-care patches

A broad-format pack for eye-care patches, shown from the printed sealing face and the formed reverse.



Thermoformed cup

Sheet mask

A pack with a broad sealing perimeter for a sheet mask. The printed face completes the single-use presentation.



Multi-compartment pack

Ketchup and mayonnaise

Ketchup and mayonnaise in two separate cavities on one pack. The central sealed area keeps the fills separate.



Multi-compartment pack

Cleanse and moisturize set

Cleanser and moisturizer in two separate cavities on one pack, combining two care steps in one set.



Multi-compartment pack

Four-sauce portion tray

Four wells in one opaque tray for tomato sauce, mayonnaise, mustard and pesto. One shared foil lid opens the pack; front and underside views.



Thermoformed cup

Ketchup portion

A single-cavity portion pack for ketchup, with product information on the printed lidding film.



Thermoformed cup

Mustard portion

A single-portion mustard pack. An edge extension provides a grip for the lidding film.



Thermoformed cup

Skincare cream

Skincare cream in a custom-formed cavity. The print area and the cavity share one integrated pack.



Thermoformed cup

Honey portion

A custom-moulded single-cavity pack for honey, presented as a portion with a film lid.



Thermoformed cup

Clay mask portion

A shallow single-cavity application for a clay mask, completed by the printed sealing face.



Thermoformed cup

Skincare cream portion

A custom-moulded pack for skincare cream. The same cavity structure is also shown with honey as an alternative fill.



Thermoformed cup

Tomato sauce

Tomato sauce in an opaque red shell-shaped cup, shown open from above and from its moulded underside.



Thermoformed cup

Face cream

Face cream in an opaque pearl-white shell-shaped cup, shown open from above and from its moulded underside.



Thermoformed cup

Pesto

Pesto in a clear leaf-shaped cup, shown open from above and through its transparent moulded underside.



Thermoformed cup

Face cream

Face cream in an opaque forest-green leaf-shaped cup, shown open from above and from its moulded underside.



Thermoformed cup

BBQ sauce

BBQ sauce in an opaque charcoal star-shaped cup, shown open from above and from its moulded underside.



Thermoformed cup

Face cream

Face cream in an opaque dusty-rose star-shaped cup, shown open from above and from its moulded underside.



Thermoformed cup

Mustard

Mustard in an opaque yellow hexagonal cup, shown open from above and from the same pack's moulded underside.



Thermoformed cup

Face cream

Face cream in a clear hexagonal cup, shown open from above and through its transparent moulded underside.



Thermoformed cup

Honey

A custom-moulded pack with a film lid for single-portion presentations.



Thermoformed cup

Clay mask

A custom-moulded pack with a film lid for single-portion presentations.



Thermoformed cup

Strawberry jam

A custom-moulded pack with a film lid for single-portion presentations.



Thermoformed cup

Cleansing gel

A custom-moulded pack with a film lid for single-portion presentations.



Thermoformed cup

Mustard sauce

A custom-moulded pack with a film lid for single-portion presentations.



Thermoformed cup

Skincare cream

A custom-moulded pack with a film lid for single-portion presentations.



Thermoformed cup

Chocolate spread

A custom-moulded pack with a film lid for single-portion presentations.




Thermoformed cup

Hair mask

A custom-moulded pack with a film lid for single-portion presentations.

From mould to filling.

Tell us about your product and how it will be used. We will assess the packaging format, fill quantity and print together.



Inside the production line

Formed cavities, dosing and the filling line.

01 Product assessment

Product properties, intended use, target fill quantity and order requirements are assessed.

02 Custom mould and sample

A custom mould is produced for the required shape, dimensions and cavity layout, followed by a sample.

03 Filling and print preparation

Filling conditions and print design are prepared according to product and packaging compatibility.

04 Production planning

Production, inspection and dispatch are planned around the approved sample and technical scope.

snappax

info@setuna.com.tr

+90 531 962 80 55

snappax.com